
THUASNE

GenuStart® ROM

fr	Genouillère ligamentaire articulée avec contrôle de flexion/extension	5
en	Hinged ligament knee brace with flexion/extension control	6
de	Knieführungsohrthese mit Beuge- und Streckbegrenzung	8
nl	Ligamentaire kniebrace met scharnier met buig-/strekcontrole	10
it	Ginocchiera per legamenti con controllo di flessione/estensione	12
es	Rodillera ligamentaria articulada con control de flexión/extensión.....	14
pt	Joelheira ligamentar articulada com controlo de flexão/extensão.....	16
da	Artikuleret knæbind til ledbånd med kontrol af fleksion/ekstension.....	18
fi	Nivelletty polven nivelsidetuki, jossa on fleksion ja ekstension rajoitusmahdollisuus.....	20
sv	Artikulerat knästöd för ligamentskada med flexions-/extensionskontroll.....	21
el	Επιγονατίδα συνδέσμων αρθρωτή με έλεγχο της κάμψης/έκτασης.....	23
cs	Kolenní ortéza s kloubem umožňující nastavení flexe/extenze	25
pl	Przegubowa orteza stawu kolanowego z kontrolą zgięcia/wyprostu	27
lv	Ceļa locītavas saišu ortoze ar eņģēm un saliekšanas/atliekšanas kontroli.....	29
lt	Lankstus raiščiu antkelis su sulenkimo / ištiesimo kontrole	31
et	Painutust-sirutust reguleerida võimaldav liigenditega põlveliigese tugikaitse.....	32
sl	Pregibna kolenska opornica za nadzorovano upogibanje in iztegovanje.....	34
sk	Kolená ortéza na podporu väzov s nastavením flexie/extenzie.....	36
hu	Ízületes térdrögzítő hajlítás-nyújtás szabályozással	38
bg	Съчленена шина за връзките на коляното с контрол на съгване/разгъване.....	40
ro	Genunchieră ligamentară articulată cu controlul flexiei/extensiei.....	42
ru	Лигментарный коленный ортез с регуляцией угла сгибания и разгибания.....	43
hr	Zglobna ortoza za ligamente koljena s kontrolom fleksije/ekstenzije.....	45
zh	帶屈曲/伸展控制關節的韌帶護膝.....	47
ar	واقى ركبّة مفصلي للاربطة مع ميزة التحكم بالتي والتمد.....	51

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال

		(a)	(b)	(c)
	XS	33 - 39 cm	30 - 33 cm	25 - 30 cm
	S	39 - 47 cm	33 - 36 cm	30 - 36 cm
	M	47 - 53 cm	36 - 38 cm	36 - 41 cm
	L	53 - 60 cm	38 - 43 cm	41 - 46 cm
	XL	60 - 67 cm	43 - 48 cm	46 - 51 cm
	XXL	67 - 75 cm	48 - 53 cm	51 - 56 cm
	XXXL	75 - 81 cm	53 - 58 cm	56 - 61 cm

fr	(a) Circonférence de la cuisse (15 cm au-dessus du genou) (b) Circonférence du genou (c) Circonférence du mollet (15 cm au-dessous du genou)
en	(a) Thigh circumference (15 cm above the knee) (b) Knee circumference (c) Calf circumference (15 cm below the knee)
de	(a) Oberschenkelumfang (15 cm über dem Knie) (b) Knieumfang (c) Wadenumfang (15 cm unter dem Knie)
nl	(a) Omvang van de dij (15 cm boven de knie) (b) Omvang van de knie (c) Omvang van de kuit (15 cm onder de knie)
it	(a) Circonferenza della coscia (15 cm al di sopra del ginocchio) (b) Circonferenza del ginocchio (c) Circonferenza del polpaccio (15 cm al di sotto del ginocchio)
es	(a) Circunferencia del muslo (15 cm por encima de la rodilla) (b) Circunferencia de la rodilla (c) Circunferencia de la pantorrilla (15 cm por debajo de la rodilla)
pt	(a) Circunferência da coxa (15 cm acima do joelho) (b) Circunferência do joelho (c) Circunferência da barriga da perna (15 cm abaixo do joelho)
da	(a) Lårets omkreds (15 cm over knæet) (b) Knæets omkreds (c) Læggets omkreds (15 cm under knæet)
fi	(a) Reiden ympärys (15 cm polven yläpuolelta) (b) Polven ympärys (c) Pohkeen ympärys (15 cm polven alapuolelta)
sv	(a) Omkrets runt låret (15 cm ovanför knäet) (b) Omkrets runt knä (c) Omkrets runt vad (15 cm under knäet)
el	(a) Περίμετρος του μηρού (15 εκ. πάνω από το γόνατο) (b) Περίμετρος του γόνατος (c) Περίμετρος της γάμπας (15 cm κάτω από το γόνατο)
cs	(a) Obvod stehna (15 cm nad kolenem) (b) Obvod kolena (c) Obvod lýtky (15 cm pod kolenem)
pl	(a) Obwód uda (15 cm powyżej kolana) (b) Obwód kolana (c) Obwód łydki (15 cm poniżej kolana)
lv	(a) Augšstilba apkārtmērs (15 cm virs ceļgala) (b) Ceļgala apkārtmērs (c) Apakšstilba apkārtmērs (15 cm zem ceļgala)
lt	(a) Šlaunies apimtis (15 cm virš kelio) (b) Kelio apimtis (c) Blauzdos apimtis (15 cm žemiau kelio)
et	(a) Reie ümbermõõt (15 cm põlvest kõrgemal) (b) Põlve ümbermõõt (c) Sääremarja ümbermõõt (15 cm põlvest allpool)
sl	(a) Obseg stegna (15 cm nad kolenom) (b) Obseg kolena (c) Obseg goleni (15 cm pod kolenom)

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Conserver cette notice.

en

HINGED LIGAMENT KNEE BRACE WITH FLEXION/EXTENSION CONTROL

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Bilateral product.

The product is available in 4 forms: 2 versions (open and closed), each available in 2 lengths (short 30 cm and long 40 cm) to meet the needs of patients (example of open long version in figure A and closed short version in figure B).

Each version is available in 7 sizes.

The device is composed of:

- a fabric part consisting of a 3D knit at the front and a thin, elastic knit at the back,
- 2 rigid hinged uprights providing knee stability (A),
- the TM5+ hinge mimicking the knee's natural movement (B), with range of motion control, protected by a cover,
- 2 anterior half-straps (②, ④) and 2 posterior half-straps (①, ③) (including a suspension strap ①) for the short version,
- 4 anterior half-straps (②, ④, ⑤, ⑦) and 3 posterior half-straps (①, ③, ⑥) (including a suspension strap ①) for the long version,
- a patellar insert to enable correct positioning of the product on the leg (C).

The closed version has finger loops for easier fitting.

The flexion/extension adjustment system is designed to be used without tools.

The design of the upright enables limitation of extension to 0° by default.

A bag containing 4 metal pins (figure C), required to adjust the hinge of the knee brace, is included in the packaging.

Composition

Textile components: polyamide - polyester - synthetic rubber - elastane - polyvinyl chloride.

Rigid components: polyoxymethylene - aluminium - stainless steel.

Properties/Mode of action

Stabilization of knee joint ligaments provided by the rigid hinged side uprights.

The brace is maintained on the leg thanks to the Suspension Strap ①.

- Extension can be adjusted to 10° and 20°.
- Flexion can be adjusted to 0°, 45°, 60° and 90°.

Indications

Conservative treatment of knee ligament injuries and/or ruptures (cruciate and/or lateral ligaments).

Post-operative rehabilitation.

Joint instability/laxity.

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

History of venous or lymphatic disorders.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

If the measurements taken do not correspond to the same size, favor the size corresponding to measurement a.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

Do not wear the product in a medical imaging machine.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.

A. Positioning of the product on the leg

For the open version:

1. Position the open product on the back of the leg, placing the center of the hinges level with the patella.
2. Cover the leg with the knee brace by fastening the touch fasteners at the top and bottom.

For the closed version:

Slip the knee brace onto the leg, making sure the patella is properly centered in the patellar comfort zone.

B. Fastening the straps

1. First tighten and fasten the suspension strap (strap ①) located above the calf.
2. Then successively tighten and fasten straps ②, ③ and ④.
3. For the long version, then tighten and fasten straps ⑤, ⑥ and ⑦.

Tightening the front/back straps adjusts the position of the uprights on the leg to optimize comfort.

If you feel the need to tighten the straps during an activity, make sure you tighten them following the steps indicated in the fitting instructions.

C. Flexion/extension control:

The adjustment of the flexion/extension must be decided and performed by the healthcare professional, not by the patient.

Extension can be limited to 10° and 20°.

Flexion can be limited to 0°, 45°, 60° and 90°.

The product's default position is limitation of flexion to 90°.

To change this setting, please follow the following instructions, **to be repeated with the same limitation on both hinges.**

1. Open the hinge cover. To remove hinge cover, push the fastener from underneath the hinge and peel back upon release. (figure D).
2. Remove the metal pin initially inserted (if necessary).
3. Insert the metal pins (figure C) according to the desired limitation angles.
The limitation angles are indicated on the sticker on the center of the hinge (figure E).
4. Close the cover again: hook the anchor under the front of the hinge and snap the fastener back in place (figure F).

The hinge cap can seem difficult to close, do not hesitate to snap it to secure.

Then perform a few flexion/extension movements to make sure limitation is properly locked at the desired angle.

To make the knee brace more comfortable and improve its hold, the uprights can be shaped.

- Apply a slight pressure in the desired direction.
- Repeat the process on the other upright if needed.

If subsequent shaping of the uprights is needed, consult a healthcare professional.

Care

Close the self-fastening tabs before washing. Hand wash. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2019103 (2022-12)